

## Green Max “Regular” Fiber Optic Reusable Blades

Product Model : **LSC22**



**Category: Medical Equipment, Laryngoscope Set & Spares,**

### Summary :

Utilizing advanced LED light technology, **GWS Green Max Fiber Optic laryngoscopes** deliver superior illumination for safer intubations. Reusable Fiber Optic blades are manufactured with 100% high quality stainless steel to ensure durability and long life. Each blade is MATT finishes to reduce glare and enhance ease of intubation for optimal patient safety.

### Specification :

- Enlarged 3.6mm diameter fiber optic bundle.
- Engineered with a removable fiber optic bundle for easy replacement of damaged bundles.
- Autoclavable Light Guide, tested to withstand over 4000 cycles.
- Matt Finish to avoid any reflection during intubation.
- Compatible to all ISO Handles 7376-3
- Available in Macintosh & Miller Blades.