


## Electric Syringe Needle Destroyer-ABS Body

Product Model : **ESYD12**



**Category: Syringe Needle Destroyer, Medical Equipment,**


### Specification :

- GWS Brand Electric Needle & Syringe Destroyer
- Electrical instant needle burner and replaceable cutter to cut syringe nozzles after use.
- 230V AC +- 10% 50 hz., single phase (alternative voltages available), protection 2 amperes 3/4
- Weight (approx.): - 1.8kgs
- Dimension(with packing):- 23 cms x 20 cms x 16 cms

### Optional:-

- Spare Cutter for syringe destroyer(ESYD12.01)

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com