

## Endotracheal Tube Plain

Product Model : **MD11**



**Category: Medical Equipment, Medical Disposables, Anaesthesia,**

### Specification :

- Designed to be inserted into the trachea through the oral/nasal cavity to secure the airway and maintain the gas exchange in the lungs.
- Made from non-toxic, non-irritant medical grade PVC
- Radio opaque line provided throughout the length for X-ray visualization
- Universal 15mm Connector at Proximal End
- The tip is suitable for both Nasal & Oral intubation
- Plain Sizes (mm) : 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10mm

### Endotracheal Tube Plain - Sizes/ Codes

Size	2.0mm	2.5mm	3.0mm	4.0mm	4.5mm	5.0mm	5.5mm	6.0mm
Code	MD11.2	MD11.2.5	MD11.3	MD11.4	MD11.4.5	MD11.5	MD11.5.5	MD11.6
Size	6.5mm	7.0mm	7.5mm	8.0mm	8.5mm	9.0mm	9.5mm	10mm

Code	MD11.6.5	MD11.7	MD11.7.5	MD11.8	MD11.8.5	MD11.9	MD11.9.5	MD10
------	----------	--------	----------	--------	----------	--------	----------	------

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)