

## Guedel Airways

Product Model : **Default Model No**



**Category: Medical Equipment, Medical Disposables, Anaesthesia,**

### Summary :

Guedel Airways are essential medical devices designed to maintain an open airway in patients who are unconscious or semi-conscious. Made from high-quality, durable materials, these airways feature a curved design that helps prevent airway obstruction by lifting the tongue away from the throat. Available in various sizes to accommodate different patient needs, Guedel Airways are easy to insert and remove, making them a crucial tool in emergency care, anesthesia, and critical care settings. Their reliability and effectiveness in ensuring clear airways make them indispensable in hospitals, clinics, and other medical environments.

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com