

## Siliconised Endotracheal Tube Plain & Cuffed

Product Model : **MD13** | **MD14**



**Category: Medical Equipment, Medical Disposables, Anaesthesia,**

### Specification :


Plain | MD13

Cuffed | MD14

- Used to maintain a patent airway and mechanically ventilate patients
- Made from a special non-toxic Silicone and PVC blend
- PHTHALATE FREE - minimizing patient exposure to carcinogenic Phthalates
- Extra Smooth, non-adherent inside/outside surface for ease in intubation, extubation and suctioning
- Kink resistant, thin-walled, thermo-sensitive tube, softens at body temperature
- Smooth polished murphy eye, distal end and tapered balloon end for atraumatic intubation and extubation
- Kink resistant inflation tube to ensure patient safety during cuff inflation and deflation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use
- Cuffed Sizes (mm) : 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10
- Plain Sizes (mm) : 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10

*Legacy of Excellence in Health Care...*

---

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmmed.com](mailto:info@gwsmmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmmed.com](http://www.gwsmmed.com)