

## Levin's Tube

Product Model : **Default Model No**



**Category: Medical Equipment, Medical Disposables, Gastro-Entrology,**

### Summary :

Levin Tube is a versatile medical device designed for gastrointestinal drainage and suction. Constructed from high-quality, flexible materials, it ensures smooth and reliable performance. The Levin Tube features multiple drainage eyes along its length to effectively remove gastric contents or fluids. Its design includes a tapered tip for easy insertion and minimal discomfort. Suitable for use in hospitals, clinics, and medical facilities, this tube is essential for managing gastrointestinal conditions and performing various diagnostic and therapeutic procedures.

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com