

## Examination Table with Height Adjustment, Hydraulic

Product Model : **HF145**



**Category: hospital furniture, Medical Equipment, Examination Couch & Tables,**

### Specification :

- Two Section Top
- MS Framework, pre-treated & epoxy powder coated.
- Adjustable Height by Hydraulic Pump
- Weight Load Capacity 180-200Kgs
- Height Adjustment(Lo-Hi):- 51cms to 82cms
- Product Dimensions:- 72? x 24? x 20?
- Also Available in Electric:- **EHF145**

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)