

Sterilization Tray with Silicone Mat, Stainless Steel

Product Model : **OST11**




Category: Stainless Steel Hospital Hollowware, Medical Equipment, Sterilization Trays,

Specification :

- Sterilization Tray, used for ophthalmic instruments sterilizing.
- Made with Stainless Steel Material.
- Silicone Mat inside.
- Press button lock for opening and closing

| Item Code | Type | Size(cms) LxWxH | Net Weight |
|-----------|--------|-----------------|------------|
| OST11.01 | Small | 20 x 6.5 x 2.5 | 230 Grams |
| OST11.02 | Medium | 21 x 12.5 x 2.5 | 400 Grams |
| OST11.03 | Large | on request | on request |

Legacy of Excellence in Health Care..

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com