

## Sterilization Tray with Silicone Mat, Stainless Steel

Product Model : **OST11**




**Category: Stainless Steel Hospital Hollowware, Medical Equipment, Sterilization Trays,**

### Specification :

- Sterilization Tray, used for ophthalmic instruments sterilizing.
- Made with Stainless Steel Material.
- Silicone Mat inside.
- Press button lock for opening and closing

Item Code	Type	Size(cms) LxWxH	Net Weight
OST11.01	Small	20 x 6.5 x 2.5	230 Grams
OST11.02	Medium	21 x 12.5 x 2.5	400 Grams
OST11.03	Large	on request	on request

*Legacy of Excellence in Health Care..*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com