

## **Black Resuscitator Autoclavable (Adult)**

Product Model : **AN-106**




**Category: Anaesthesia Products & Equipment, Medical Equipment, Artificial Resuscitators(Ambu Bag) Black Rubber,**

### **Specification :**

- Double Layer Autoclavable Bellow.
- Black Rubber Super Face Mask Size: 4 no.
- Autoclavable & Reusable
- Zip Pouch Packing
- Color-Black
- Packing: Zip Pouch Packing

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwsurgicals@gmail.com](mailto:gwsurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)