

## Black Resuscitator Autoclavable (Child)

Product Model : **AN-107**



**Category: Medical Equipment, Anaesthesia Products & Equipment, Artificial Resuscitators(Ambu Bag) Black Rubber,**

### Specification :

- Double Layer Autoclavable Bellow.
- Black Rubber Super Face Mask Size: 2 no.
- Autoclavable & Reusable
- Zip Pouch Packing
- Color-Black
- Packing: Zip Pouch Packing

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com