

## RESPO-Infant Resuscitator

Product Model : **RESP01**



**Category: Baby Care Equipment, Medical Equipment, Neonatal Resuscitation System,**


### Specification :

- Used on delivery room, baby wards, NICU that the baby's maximum weight up to 10KG
- Inside an alarm system of differential pressure. When alarming, the safety valve starts automatically and produces the mixture gas continuously to ensure the safety of patients
- The range of flow adjustment is 5LPM~15LPM.
- The MAX-P and PIP are safe, controlled and accurate, PEEP is constant and accurate.

### Technical Specifications:-

- Pressure gauge: -2~10kPa (-20~100cmH<sub>2</sub>O).
- MAX-P: @0~15LPM 6kPa (60cmH<sub>2</sub>O).
- PIP: @5~15LPM 0.2~5.7kPa (2~57cmH<sub>2</sub>O).
- PEEP: @5~15LPM 0.2~5.7kPa (2~57cmH<sub>2</sub>O).

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com