

RESPO(PRO)-Infant Resuscitator

Product Model : **RESP02**



Category: Baby Care Equipment, Medical Equipment, Neonatal Resuscitation System,


Specification :

- Used in delivery room, baby wards, NICU that the baby's maximum weight up to 10KG
- Inside an alarm system of differential pressure. When alarming, the safety valve starts automatically and produces the mixture gas continuously to ensure the safety of patients
- Adjust the Oxygen concentration from 21% to 100% optionally
- The range of flow adjustment is 5LPM~15LPM.
- The MAX-P and PIP are safe, controlled and accurate, the PEEP is constant and accurate.

Technical Specifications:-

- Pressure gauge: -2~10kPa (-20~100cmH₂O).
- MAX-P: @0~15LPM ?6kPa (60cmH₂O).
- PIP: @5~15LPM @ 0.5~1.5kPa (5~15cmH₂O).

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwsurgicals@gmail.com  +91-124-4915548  www.gwsmed.com