

## Ophthalmic Operation Table, Hydraulic

Product Model : **OPH13**



**Category: Operation Table, Medical Equipment, Ophthalmic Operation Table,**

### **Specification :**

- Ophthalmic OT Table has a maximum height of 920 mm and minimum height of 690 mm.
- Ideal for standing table or sitting working position.
- It is an ideal operating table with the advantages of strong applicability, easy to operate.
- Adjusts conveniently and steady performance.
- GWS-OPH13 its ergonomic design ensures both surgeon and patient comfort.
- Patient comfort is assured with cushions which are easy to clean.
- The adjustable cushioned head-rest ensures doctor comfort and resistant during operations.

### **Technical Specification:-**


- Top Dimension: 1900 mm (L) x 600mm (W)
- Width of Table Top 520mm
- Minimum Height: 690 mm
- Maximum Height: 920 mm
- Stroke length: 220 mm
- Patient Weight Capacity: 170 kgs.

- Back rest section: 20? (head side)
- Leg Section: 30?

Optional:- **Electric Model(OPH13-E)**

- Power Requirement: 230V AC /50Hz, 110V AC /60 Hz
- Fuses: 2A Slow Blow
- Power in VA: 150VA

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)