

## **Baby Bassinet with utility box**

Product Model : **HF-210**




**Category: hospital furniture, Medical Equipment, Baby Crib & Cradles,**

### **Specification :**

- Frame made up of ERW tubular pipe of 2.5cms. dia. Mounted on four 5cms dia. castors.
- Removable CRIB of poly glass titable to trendelenburg and R-trend positions.
- Laminated bottom shelf with a spacious cabinet for storage.
- Provision for monitor stand and drip rod available.
- Pre-treated and epoxy powder coated.

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwsurgicals@gmail.com](mailto:gwsurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)