

## Stretcher on Trolley

Product Model : **HF-133**



**Category: hospital furniture, Medical Equipment, Stretchers,**

### Specification :

- Removable with fixed height.
- Frame is made of CRCA tube.
- Trolley mounted over 150mm castors.
- IV/Saline stand Provision.
- Dimension(approx.): 2100 L x 560 W x 810 H(mm)

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwsurgicals@gmail.com](mailto:gwsurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)