

Bed Side Screen (4 fold)

Product Model : **HF-200**



Category: hospital furniture, Medical Equipment, Bed Side Screen,

Specification :

- Made of CRCA pipes on 5cms castors provided with-out curtain / springs. Powder coated.
- Dimension: 66" H x 96" L
- * Bed Side Screen Three Fold (Code No: HF-200.01)

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwsurgicals@gmail.com  +91-124-4915548  www.gwsmed.com