

## Obstetrics / Labour Table

Product Model : **HF-136**



**Category: hospital furniture, Medical Equipment, Obstetric Labour Tables,**

### Specification :

- Three section top made up of SS.
- Leg section sliding completely under the main section.
- Base frame made up of CRCA tubes.
- Trendlenburg's position achieved by pneumatic pump.
- Backrest adjustable with ratchet mechanism.
- Lithotomy crutches / Lithotomy rods with straps.
- Pre-treated & epoxy powder coated.
- Dimension(approx.): 860 L x 680 W x 760 H (mm)
- Optional:- Mattress (HF-136.01)

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com