

## Overbed Table (Fix Height)

Product Model : **HF-160**



**Category: hospital furniture, Medical Equipment, Overbed Tables,**

### Specification :

- Tubular framework mounted on 5cm castors.
- Fixed wooden laminated top.
- Pre-treated & epoxy powder coated.
- Dimension(apprx.): 1120 L x 400 W x 990 H (mm)

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com