

## Double Bowl Stand

Product Model : **HF-189**




**Category: hospital furniture, Medical Equipment, Bucket, Basins & stands,**

### Specification :

- Made with MS tubular frame
- Two bowl fitted in ring of 35cm
- Supplied with two litre stainless steel bowl
- Mounted on 5cm dia. castors
- Pre-treated and epoxy powder coated finish

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com