

## HI-VAC, MS Frame

Product Model : **SU-101**




**Category: Suction Machine, Medical Equipment, Electric Suction Units,**

**Specification :**

Housing	MS Powder coated on four castors
Max. Vacuum/LPM	-700mm Hg, 50 Ltr. / min.
Pump Type	Double Rotary Vane Pump
Jars	2 x 2 Ltr. PC Jar
Filter	Trap Filter
Tubing	10 mm ID x 2 mtr. (PVC)
Vacuum Gauge	10 cm dia. 0-760 mm Hg.
Power Supply	220 V AC, 50 / 60 Hz , 380 watt
Noise Level	50 dB A $\pm$ 3
Dimensions	49 x 36 x 70 cm
Net Weight	37 Kg
Overflow Safety	Mechanical

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmmed.com](mailto:info@gwsmmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmmed.com](http://www.gwsmmed.com)