

## Dual Power Suction Unit (Electric cum Manual)

Product Model : **SU-21**




**Category: Suction Machine, Medical Equipment, Electric cum Manual Suction Units,**

### Specification :

- **Housing:** MS Powder Coated, SS Top on two wheels
- **Capacity (Electrical):** -600 mm Hg  $\pm$  10 , 45 Ltrs/ min
- **Capacity (Pedal):** -600 mm Hg  $\pm$  10 , 200ml/ Push Stroke
- **Pump Type:** Electric – Diaphragm / Pedal – Piston
- **Jars:** 2 x 2 Ltrs Polycarbonate Jar
- **Filter:** Bacterial Filter Autoclavable/Reusable
- **Tubing:** 10mm ID x 2 mtr. (PVC)
- **Vacuum Gauge:** 6.25 cm dia 0-760 mm Hg calibration.
- **Power:** 220 V AC, 50 / 60 Hz, 90Watt (110V on request)
- **Noise Level:** 55 dB A  $\pm$  3
- **Dimension & Weight:** 42 x 39 x 64 cms, 15.5 Kg
- **Overflow Safety:** Dual Mechanical

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com