

Electric Cum Foot operated Suction

Product Model : **SU-23**



Vacuum Pump



600ml PC Jar




Oil Free

Category: Suction Machine, Medical Equipment, Electric cum Manual Suction Units,

Specification :

- Housing: Powder Coated anti corrosive MS
- Jars: 600ml Poly-carbonate
- Overflow protection: Mechanical type
- Tubing: Non-collapsible suction tubing.
- Vacuum gauge: 2.0 inch, 0-760mm Hg
- Filter: Autoclavable bacterial filter.
- Pump type: Oil free piston & manual pump
- Capacity: -580mm/ Hg +/- 10 at 12~15 LPM
- Power: 220/ 230V AC / 12V DC
- Backup: 120 minutes

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwsurgicals@gmail.com  +91-124-4915548  www.gwsmed.com