

## **Zono Refraction Chair**

Product Model : **OPH-02**



**Category: Ophthalmic Instruments, Refraction Chair,**

### **Specification :**

#### **Features:**

- Smooth-rotating, well-balanced Sliding Table Top with ample space for two instruments
- Automatic, under-the-table connection for Slit Lamp
- Touch-activated Safety Switch to stop inadvertent chair movements
- Integrated, centralized Control Console with touchstone keys
- Strong Pole to support Lamp, Phoropter, Projector & Near Vision Holder
- Convenient Instrument Holder Tray
- Dust-free, enclosed Trial Lens Rack with straight, unobtrusive, in & out sliding movement.

#### **CONTROL PANEL:**

- Elegant, single-piece Console panel for centralized unit control
- Multifaceted controls for chair movements
- Built-in connectivity for Indirect Ophthalmoscope, Near Vision Holder & Projector.

## CHAIR SPECIFICATION

- Fully flexible, rotating recliner chair
- Silent, CE-marked, German DC motor-driven, up & down movements and inclination, and Zero Position Memory
- Adjustable footrest, which can be fully raised to sleeping position
- Foot Switch for fluid up / down movement control
- Comfortable polyurethane-filled Seat
- Full movement control from backrest panel.

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com