

## 3.5mm Small Cannulated Cancellous Screw, Fully Threaded

Product Model : **2200**



**Category: Cannulated Cancellous Screw, Orthopaedic Implants,**


**Specification :**

Technical Specification		
	Thread Dia	3.5mm
	Drill Bit for Threaded Hole	2.5mm
	Core Dia.	1.8mm
	Hexagonal Socket Width	2.5mm
	Head Dia.	6.0mm
	Thread Length	Fully
	Lengths	10mm to 60mm

Codes (S.S)	Length(mm)	Codes (Titanium)
2200.35.10	10 mm	Ti-2200.35.10
2200.35.12	12 mm	Ti-2200.35.12

2200.35.14	14 mm	Ti-2200.35.14
2200.35.16	16 mm	Ti-2200.35.16
2200.35.18	18 mm	Ti-2200.35.18
2200.35.20	20 mm	Ti-2200.35.20
2200.35.22	22 mm	Ti-2200.35.22
2200.35.24	24 mm	Ti-2200.35.24
2200.35.26	26 mm	Ti-2200.35.26
2200.35.28	28 mm	Ti-2200.35.28
2200.35.30	30 mm	Ti-2200.35.30
2200.35.32	32 mm	Ti-2200.35.32
2200.35.34	34 mm	Ti-2200.35.34
2200.35.36	36 mm	Ti-2200.35.36
2200.35.38	38 mm	Ti-2200.35.38
2200.35.40	40 mm	Ti-2200.35.40
2200.35.42	42 mm	Ti-2200.35.42
2200.35.44	44 mm	Ti-2200.35.44

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com 
  +91-124-4915548 
  www.gwsmed.com