

3.5mm Small Cannulated Cancellous Screw, Fully Threaded

Product Model : **2200**



Category: Cannulated Cancellous Screw, Orthopaedic Implants,


Specification :

Technical Specification		
Thread Dia		3.5mm
Drill Bit for Threaded Hole		2.5mm
Core Dia.		1.8mm
Hexagonal Socket Width		2.5mm
Head Dia.		6.0mm
Thread Length		Fully
Lengths		10mm to 60mm

Codes (S.S)	Length(mm)	Codes (Titanium)
2200.35.10	10 mm	Ti-2200.35.10
2200.35.12	12 mm	Ti-2200.35.12

2200.35.14	14 mm	Ti-2200.35.14
2200.35.16	16 mm	Ti-2200.35.16
2200.35.18	18 mm	Ti-2200.35.18
2200.35.20	20 mm	Ti-2200.35.20
2200.35.22	22 mm	Ti-2200.35.22
2200.35.24	24 mm	Ti-2200.35.24
2200.35.26	26 mm	Ti-2200.35.26
2200.35.28	28 mm	Ti-2200.35.28
2200.35.30	30 mm	Ti-2200.35.30
2200.35.32	32 mm	Ti-2200.35.32
2200.35.34	34 mm	Ti-2200.35.34
2200.35.36	36 mm	Ti-2200.35.36
2200.35.38	38 mm	Ti-2200.35.38
2200.35.40	40 mm	Ti-2200.35.40
2200.35.42	42 mm	Ti-2200.35.42
2200.35.44	44 mm	Ti-2200.35.44

Legacy of Excellence in Health Care..

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com
  +91-124-4915548
  www.gwsmed.com