

Perforated Instrument Tray without Cover

Product Model : **HH25**





Category: Stainless Steel Hospital Hollowware, Medical Equipment,

Specification :

- Made with High Quality Stainless Steel.
- Standard Height 2" inches

| Item Code | Sizes |
|-----------|-----------|
| HH25.01 | 18" X 12" |
| HH25.02 | 15" X 12" |
| HH25.03 | 12" X 10" |
| HH25.04 | 12" X 8" |
| HH25.05 | 11" X 7" |

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com
  +91-124-4915548
  www.gwsmed.com