

## Perforated Sterilization Tray with Silicone Mat

Product Model : **HH26**




**Category: Stainless Steel Hospital Hollowware, Medical Equipment,**

### Specification :

- Made from high quality 308/18/8 Stainless Steel

Item Code	Sizes (inches)
HH26.01	8" X 4" X 1"
HH26.02	9" X 8" X 1"
HH26.03	12" X 8.3" / 4 X 2"
HH26.04	14" X 8" X 1"
HH26.05	16" X 9.5" X 1"

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com 
  +91-124-4915548 
  www.gwsmed.com