

6.5mm Cannulated Cancellous Screw, Hexagonal Socket-16mm Threaded

Product Model : **2209**




Category: Cannulated Cancellous Screw, Orthopaedic Implants,

Specification :

- Thread Length : 16mm
- Drill Bit for Threaded Hole : 3.2mm
- Shaft Dia. : 4.5mm
- Drill Bit for Shaft : 4.5mm
- Length:- 50-130mm

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com