

## 6.5mm Cannulated Cancellous Screw, Hexagonal Socket-16mm Threaded

Product Model : **2209**




**Category: Cannulated Cancellous Screw, Orthopaedic Implants,**

### Specification :

- Thread Length : 16mm
- Drill Bit for Threaded Hole : 3.2mm
- Shaft Dia. : 4.5mm
- Drill Bit for Shaft : 4.5mm
- Length:- 50-130mm

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)