

6.5mm Cannulated Cancellous Screw, Hexagonal Socket-Fully Threaded

Product Model : **2211**




Category: Cannulated Cancellous Screw, Orthopaedic Implants,

Specification :

- Available in S.S & Titanium Material.

Thread Dia	6.5 mm
Head Dia	8.0 mm
Core Dia	4.5 mm
Lengths	50-130 mm
Guide wire used	1.8 mm
Cannulated Drill Bit	4.5 mm

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com