

## 6.5mm Cannulated Cancellous Screw, Hexagonal Socket-Fully Threaded

Product Model : **2211**



**Category: Cannulated Cancellous Screw, Orthopaedic Implants,**

### Specification :

- Available in S.S & Titanium Material.

<b>Thread Dia</b>	6.5 mm
<b>Head Dia</b>	8.0 mm
<b>Core Dia</b>	4.5 mm
<b>Lengths</b>	50-130 mm
<b>Guide wire used</b>	1.8 mm
<b>Cannulated Drill Bit</b>	4.5 mm

*Legacy of Excellence in Health Care...*