

3.5mm Small Cannulated Cancellous Screw, Partially Threaded

Product Model : **2201**



Category: Cannulated Cancellous Screw, Orthopaedic Implants,

Specification :

- Available in S.S & Titanium Material.

Technical Specification		
	Thread Dia	3.5mm
	Drill Bit for Threaded Hole	2.5mm
	Core Dia.	1.9mm
	Hexagonal Socket Width	2.5mm
	Head Dia.	6.0mm
	Thread Length	Partially

Legacy of Excellence in Health Care...