

3.5mm Small Cannulated Cancellous Screw, Partially Threaded

Product Model : **2201**




Category: Cannulated Cancellous Screw, Orthopaedic Implants,

Specification :

- Available in S.S & Titanium Material.

Technical Specification		
	Thread Dia	3.5mm
	Drill Bit for Threaded Hole	2.5mm
	Core Dia.	1.9mm
	Hexagonal Socket Width	2.5mm
	Head Dia.	6.0mm
	Thread Length	Partially

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com