

## **4.0mm Cannulated Cancellous Screw, Hexagonal Socket-Partially Threaded**

Product Model : **2202**



**Category: Cannulated Cancellous Screw, Orthopaedic Implants,**

**Specification :**

- Available in S.S & Titanium Material

Technical Specification		
	Thread Dia	4.0mm
	Drill Bit for Threaded Hole	2.5mm
	Shaft Dia	2.4mm
	Core Dia.	1.9mm
	Hexagonal Socket Width	2.5mm
	Head Dia.	6.0mm
	Thread Length	Partially
	Lengths	10mm to 80mm

*Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com