

## **4.0mm Cannulated Cancellous Screw, Hexagonal Socket-Partially Threaded**

Product Model : **2202**




**Category: Cannulated Cancellous Screw, Orthopaedic Implants,**

### **Specification :**

- Available in S.S & Titanium Material

Technical Specification	
Thread Dia	4.0mm
Drill Bit for Threaded Hole	2.5mm
Shaft Dia	2.4mm
Core Dia.	1.9mm
Hexagonal Socket Width	2.5mm
Head Dia.	6.0mm
Thread Length	Partially
Lengths	10mm to 80mm

*Legacy of Excellence in Health Care..*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com 
  +91-124-4915548 
  www.gwsmed.com