

4.0mm Cannulated Cancellous Screw, Hexagonal Socket-Partially Threaded

Product Model : **2202**




Category: Cannulated Cancellous Screw, Orthopaedic Implants,

Specification :

- Available in S.S & Titanium Material

| Technical Specification | |
|-----------------------------|--------------|
| Thread Dia | 4.0mm |
| Drill Bit for Threaded Hole | 2.5mm |
| Shaft Dia | 2.4mm |
| Core Dia. | 1.9mm |
| Hexagonal Socket Width | 2.5mm |
| Head Dia. | 6.0mm |
| Thread Length | Partially |
| Lengths | 10mm to 80mm |

Legacy of Excellence in Health Care..

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
 Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com
  +91-124-4915548
  www.gwsmed.com