

## Isolation Gown

Product Model : **DISP23**



**Category: Non Woven, Medical Disposables, Medical Equipment, Covid-19,**

### Summary :


Isolation Gown is a medical protective clothing made of highly effective filtering material. They are compact in structure and can effectively block liquids and particles. At the same time, they are covered with a special membrane that is inspired by protective clothing. Consistent electricity, anti-infiltration, while resistant to a variety of chemicals, liquids, acids and alkalis, and good comfort. Our products are based on the nature of the application itself. The straps of the product are elastically designed to make them more secure and the elastic band makes us flexible and easy to work with. Our products are widely used in clean workshops in hospitals, electronics, pharmacies, food, laboratories, and other industries.

### Specification :

- Material: PP or PP+PE
- GSM: 30/35/40/60
- Color Option: Blue/ Yellow
- Tie at Neck or Velcro at neck
- Round neck or V-neck
- Elastic or knitting cuff
- Ties at back of waist
- Side Tie Optional

*Legacy of Excellence in Health Care..*

---

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)