


Portable Anaesthesia Machine

Product Model : **AN-164**



Category: Anaesthesia Products & Equipment, Anaesthesia Machine,

Specification :

- Construction: trolley made of tubular rigid steel section electrostatically powder coated or stainless steel.
- Suitable drawer
- Unitised rotamotor block for oxygen & nitrous oxide with fine adjustment controls & colour coated knob and luminous back plate.
- Oxygen flow meter calibrated from 100cc/mm to 100 ltrs. / Minute & Nitrous oxide flow meter from 200cc/min to 12 liters/min down Stream receipt of oxygen
- One gold man type vaporizer for halothene
- Non return cum emergency oxygen unit with blow off valve in circuit ensures that pressure does not build-up beyond 200cm. water column & minimises risk of back flow of gases.
- Two pre set regulators (forged quality) one each for O₂ & N₂O with high pressure tubing non inter changeable for connection with flow meter inlet
- One Maggile & brain circuit

Legacy of Excellence in Health Care...



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com