


## Portable Anaesthesia Machine

Product Model : **AN-164**




**Category: Anaesthesia Products & Equipment, Anaesthesia Machine,**

### Specification :

- Construction: trolley made of tubular rigid steel section electrostatically powder coated or stainless steel.
- Suitable drawer
- Unitised rotamotor block for oxygen & nitrous oxide with fine adjustment controls & colour coated knob and luminous back plate.
- Oxygen flow meter calibrated from 100cc/mm to 100 ltrs. / Minute & Nitrous oxide flow meter from 200cc/min to 12 liters/min down Stream receipt of oxygen
- One gold man type vaporizer for halothene
- Non return cum emergency oxygen unit with blow off valve in circuit ensures that pressure does not build-up beyond 200cm. water column & minimises risk of back flow of gases.
- Two pre set regulators (forged quality) one each for O<sub>2</sub> & N<sub>2</sub>O with high pressure tubing non inter changeable for connection with flow meter inlet
- One Maggile & brain circuit

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwsurgicals@gmail.com](mailto:gwsurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)