

Green Max “SPR” Blades - Fiber Optic (4.6mm Light Guide)

Product Model : **LSC21**



Category: Medical Equipment, Laryngoscope Set & Spares,


Summary :

GWS Fiber Optic Laryngoscopes deliver superior illumination for safer intubations. **Reusable Fiber Optic blades** are manufactured with 100% high quality Stainless Steel to ensure durability and long life. Each blade is MATT finished to reduce glare and enhance ease of intubation for optimal patient safety.

Specification :

- Blades with non-replaceable fibre optics.
- All blades are made of premium grade stainless steel in matt finish.
- Fibre optics Ø 4.6mm / 6800 single fibres minimum.
- Environment friendly, economical & cost effective.
- Can be treated in autoclave up to 134° - approx. 4000 times.
- In accordance with DIN ISO 7376

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmmed.com