

## Green Max “SPR” Blades - Fiber Optic (4.6mm Light Guide)

Product Model : **LSC21**



**Category: Medical Equipment, Laryngoscope Set & Spares,**

### Summary :


**GWS Fiber Optic Laryngoscopes** deliver superior illumination for safer intubations. **Reusable Fiber Optic blades** are manufactured with 100% high quality Stainless Steel to ensure durability and long life. Each blade is MATT finished to reduce glare and enhance ease of intubation for optimal patient safety.

### Specification :

- Blades with non-replaceable fibre optics.
- All blades are made of premium grade stainless steel in matt finish.
- Fibre optics Ø 4.6mm / 6800 single fibres minimum.
- Environment friendly, economical & cost effective.
- Can be treated in autoclave up to 134° - approx. 4000 times.
- In accordance with DIN ISO 7376

*Legacy of Excellence in Health Care...*

---

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmmed.com](mailto:info@gwsmmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmmed.com](http://www.gwsmmed.com)