


Flexitip Optic LED Balde

Product Model : **LSC13**



Category: Medical Equipment, Laryngoscope Set & Spares,

Summary :


The Forte **GWS Fliptip blade** is specially designed to improve visualization of the larynx during intubations. With a blade tip adjustable through 80degrees the Forte Fliptip blade can simplify elevation of the epiglottis for many conditions, including reduced thyromental distance, short or immobilized neck, obesity and swollen tissue of the larynx.

Forte Fliptip laryngoscope blades are ideal for EMS, Emergency Room, and any other setting with unpredictable patient airway conditions.

Specification :

- LED technology provides up to 40% more light than conventional halogen bulbs.
- MATT finished blade surface reduces glare to enhance ease of intubation.
- Compatible with all ISO 7376-1 Conventional handles.

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com