

## **Eco Frame Back Support**

Product Model : **RHA519**



**Category: Medical Equipment, Rehabilitation Supports, Back Support,**

### **Specification :**

**Features:** Made from special tempered spring steel for rigid and comfortable back support. It immobilizes the lumbar region effectively.

**Applications:** Immobilization of the lumbar region. Used to give support to the back.

Size Inches	Item Code
28"	RHA519.1
30"	RHA519.2
32"	RHA519.3
34"	RHA519.4

36"	RHA519.5
38"	RHA519.6
40"	RHA519.7
42"	RHA519.8
44"	RHA519.9

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)