

L.S (Lumbar Sacro) Belt, Primo

Product Model : **RHA512**



Category: Medical Equipment, Rehabilitation Supports, Back Support,

Specification :

Features: Elastic supports on both sides for snug fit. Front and back panels are made of soft foam for comfort. Hook and loop closures for easy application and removal. Two rigid aluminium bars and two plastic bars for reinforcement. Double stitching on the edges for durability


Application:

- Chronic low back pain & strain.
- Osteoporotic pain in the lumbosacral region.
- Degenerative Lumbar Disc Disease.
- Post surgical lumbo sacral immobilization.

| | |
|-------------------------------|-----------|
| Sacro Lumbar (L.S) Belt Primo | |
| Sizes | Item Code |

| | |
|----------|----------|
| Small | RHA512.1 |
| Medium | RHA512.2 |
| Large | RHA512.3 |
| X Large | RHA512.4 |
| XX Large | RHA512.5 |

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com