

Ambulance Collar

Product Model : **RHA529**



Category: Medical Equipment, Rehabilitation Supports, Cervical Support,

Summary :

The Ambulance Collar is an essential piece of equipment in the event of a neck injury. This cervical collar supports and stabilizes the neck, preventing further damage and helping to ensure a speedy recovery. The Ambulance Collar is made from high-quality materials and is available at a competitive price.

Specification :

- **Features:** Single Piece Flat Design. Semi Rigid Plastic, Padded With Sturdy Foam Provides Comfortable Immobilization. It Is Used For Temporary Cervical Immobilization For Accidental Victims During The Transportation & Diagnostic Studies.
- **Application:** It Is Used As Temporary Support To Neck In Accidental Conditions. It Provides Rigid Support To Limit The Neck Movement Or To Prevent Harmful Neck Movements.

| Ambulance Collar | |
|-------------------------|-----------|
| Size | Item Code |
| Small | RHA529.1 |
| Medium | RHA529.2 |
| Large | RHA529.3 |

Legacy of Excellence in Health Care...



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com