

## Cervical Collar Hard

Product Model : **RHA534**



**Category: Medical Equipment, Rehabilitation Supports, Cervical Support,**

### Summary :

Are you looking for a quality cervical collar that is both comfortable and effective? Look no further than the Cervical Collar Hard from GWS Surgicals LLP. Made of light, high-density polythene sheet, this collar features proper foam padding for all-day comfort. Plus, it comes with velcro closure for easy application. And to ensure maximum ventilation, the collar has strategically placed holes. So if you're in need of a top-quality cervical collar, look no further than the Cervical Collar Hard from GWS Surgicals LLP.

### Specification :

**Feature:** Made of light high density polythene sheet. Proper foam padding for comfortable wear. Velcro closure for easy application. Holes for ventilation.

**Application:** Cervical spine immobilization in pre & post surgical cases. Cervical spondylosis.

Size	Item Code
Small	RHA534.1
Medium	RHA534.2
Large	RHA534.3
X-Large	RHA534.4

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com