

## Cervical Collar Hard Adjustable

Product Model : **RHA528**



**Category: Medical Equipment, Rehabilitation Supports, Cervical Support,**

### Summary :

Looking for a sturdy cervical collar that can be easily adjusted? Look no further than the Cervical Collar Hard Adjustable from GWS Surgicals LLP. This hard collar is made of high-quality materials and can be easily adjusted to provide the perfect level of support and comfort. Whether you're looking for an orthopedic support for daily use or a bulk manufacturer for export, GWS Surgicals LLP has you covered.

### Specification :

- **Features:** Made of light high density polythene sheet. Proper foam padding for comfortable wear. Double layer for height adjustment. Holes for ventilation. Velcro closure for easy application
- **Applications:** Cervical spine immobilization in pre & post surgical cases.  
Cervical spondylosis.

Cervical Collar Hard Adjustable	
Size	Item Code

Small	RHA528.1
Medium	RHA528.2
Large	RHA528.3
X-Large	RHA528.4

### *Legacy of Excellence in Health Care...*



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgicals@gmail.com



+91-124-4915548



www.gwsmed.com