

## 3D Knee Cap

Product Model : **RHA509**



**Category: Medical Equipment, Rehabilitation Supports, Knee/ Calf & Ankle Support's,**

### Summary :

The 3D Knee Cap is a knee support that helps to stabilize and protect the knee joint. It is made of breathable fabric and features a three-dimensional design that provides extra support and comfort. The cap is ideal for people who suffer from chronic pain, arthritis, or other joint injuries. Available in different sizes like Small, Medium, Large, X-Large.


### Specification :

- Used in case of Osteoporosis Knee ,Knee swelling, Knee pain.
- Anatomically shaped with “reduced compression spot” ensures no patellar pressure.
- Three dimensional knitting provides anatomical shape, snug fitting and strong compression.
- Bi-layered construction. Dermophillic cotton inside ensures comfort in all weather conditions.
- Nylon outside ensures durability & excellent aesthetics and color fastness. Thick knit retains therapeutic warmth and provides excellent support.

**3D Knee Cap**

Sizes	Item Code
Small	RHA509.1
Medium	RHA509.2
Large	RHA509.3
X-Large	RHA509.4

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)