

Knee Immobilizer 14"/19"/22"

Product Model : **RHA504 | RHA505 | RHA506**



Category: Medical Equipment, Rehabilitation Supports, Knee/ Calf & Ankle Support's,

Summary :

GWS Brand! Knee Immobilizer is perfect for immobilizing the knee joint and providing support. It is made of high-quality materials and can be used for both long and short term immobilization. Our Knee Immobilizer is also available in a variety of sizes like 14" / 19" / 22" to ensure a perfect fit.

Specification :

Features: Its Contoured posterior and sturdy bilateral metal stays, immobilizes the leg at 0°-20°. It also has a removable medial, lateral and posterior stays. It has comfortable washable, skin friendly and laminated cloth material. Its patellar opening releases patellar pressure ensuring comfort. Its five straps of hook & loop fasteners help in secure fitting and better immobilization.

Application:

- Post Total Knee Replacement.
- Post Ligament Reconstruction.

Color Option: Black/ Grey/ Beige

Knee Immobilizer / Brace 14"			
Item Code	Size	Inches	Cms
RHA504.1	Small	13-14	33-36
RHA504.2	Medium	14-15	36-38
RHA504.3	Large	15-16	38-41
RHA504.4	X-Large	16-18	41-46
RHA504.5	XX-Large	18-20	46-51

Knee Immobilizer / Brace 19"			
Item Code	Size	Inches	Cms
RHA505.1	Small	13-14	33-36
RHA505.2	Medium	14-15	36-38
RHA505.3	Large	15-16	38-41
RHA505.4	X-Large	16-18	41-46
RHA505.5	XX-Large	18-20	46-51

Knee Immobilizer / Brace 22"			
Item Code	Size	Inches	Cms
RHA506.1	Small	13-14	33-36
RHA506.2	Medium	14-15	36-38
RHA506.3	Large	15-16	38-41
RHA506.4	X-Large	16-18	41-46
RHA506.5	XX-Large	18-20	46-51

Legacy of Excellence in Health Care...



Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India



info@gwsmed.com, gwssurgical@gmail.com



+91-124-4915548



www.gwsmed.com